


APPLICATION NOTE: 321



JCAHO PC.17.10

PARTS: B6, B8, C5, & C7 - TISSUE STORAGE

- B6.** Maintain continuous temperature monitoring for storage refrigerators and freezers.
- B8.** Storage equipment has functional alarms and emergency backup.
- C5.** The organization's records including storage temperatures...are retained for a minimum of ten years, or longer if required by state and/or federal laws.
- C7.** Maintain daily records to show that tissues were stored at the required temperatures.



ThermaViewer

For JCAHO PC.17.10 compliance

ThermaViewer Compliance Features:

- Monitors two refrigerators/freezers continuously (B6)
- 9 volt battery for emergency backup (B8)
- Functional alarm – phone dialer (B8)
- Downloadable temperature history (C5)
- LCD chart and PC interface for daily records (C7)

Easy to Use:

- Installs in minutes;
- No special training required;
- No maintenance;
- Stand-alone.

APPLICATION NOTE: 321



JCAHO PC.17.10 compliant
for monitoring and documenting tissue storage

